



Future State: Bridge Transfusion - Bridge Breast Milk - Available Bottles

Cerner Workflow ID: 6974 (v. 10.0) Client Workflow ID: 532

Last updated by Reiser Sanders, Feb 06, 2024 12:24pm (UTC -4 hours)

Workflow Details:

Workflow Name: Bridge Transfusion - Bridge Breast Milk - Available Bottles

Workflow State: Future State

Workstream: Other

Venue: Acute Care

Client Owner:

Cerner Owner:

Standard: Yes

Related Workflow(s):

Tags:

Workflow Summary:

Service Line:

Related Solution(s): Bridge Breast Milk Management

Project Name: Niagara Health System:OPT-0297674:NIAG_CD Niagara HIS RFP

TestBuilder Script(s):

Cerner Workflow ID: 6974 (v. 10.0)

Client Workflow ID: 532

Workflow Notes:

Introduced By: WS 3.2

Validated By: WS 5.2

Swim Lane:

Role(s): Nurse (RN)

Department(s):

Security Position(s): Nurse - Neonatology

Women's Health - Nurse

Women's Health - Clerk

Women's Health - Lactation Consultant

Start/Stop [3]

Description: Access Cerner Bridge

Decision [22]

Description: Do you want to view one specific bottle's information?

Work Step [68]

Description: Select the desired bottle or search for the bottle by typing the bottle number into the search bar & clicking the magnifying glass

Future State: Bridge Transfusion - Bridge Breast Milk - Available Bottles

Cerner Workflow ID: 6974 (v. 10.0) Client Workflow ID: 532

Last updated by Reiser Sanders, Feb 06, 2024 12:24pm (UTC -4 hours)

Work Step [39]

Description: Click on the bottle number to view fortifier information and division information

Work Step [36]

Description: Select Back to return to the Browse Bottle History screen

Work Step [4]

Description: Select Milk Mgmt

Decision [23]

Description: Do you want to edit a bottle's information?

Work Step [69]

Description: Select the desired bottle or search for the bottle by typing the bottle number into the search bar & clicking the magnifying glass

Work Step [38]

Description: Click on the pencil icon to edit the collection date/time, storage state, and/or dispose information.

Work Step [35]

Description: Click Confirm to save changes

Start/Stop [42]

Description: Select Back to return to Milk Mgmt

Work Step [5]

Description: Select Available Patient Bottles

Decision [24]

Description: Do you want to reprint a bottle label?

Work Step [30]

Description: Select the desired bottle or search for the bottle by typing the bottle number into the search bar & clicking the magnifying glass

Work Step [37]

Description: Click the Print hyperlink

Future State: Bridge Transfusion - Bridge Breast Milk - Available Bottles

Cerner Workflow ID: 6974 (v. 10.0) Client Workflow ID: 532

Last updated by Reiser Sanders, Feb 06, 2024 12:24pm (UTC -4 hours)

Work Step [34]

Description: Select the desired printer by scanning its barcode or selecting it from the list

Work Step [6]

Description: Select the desired POC and Patient Status

Decision [26]

Description: Do you want to discard a bottle?

Work Step [119]

Description: Select the desired bottle or search for the bottle by typing the bottle number into the search bar & clicking the magnifying glass

Work Step [40]

Description: Click on the dispose checkbox

Work Step [15]

Description: View the available bottle counts for the search parameters entered in the previous step

Decision [25]

Description: Do you want to view multiple bottles by sorting them?

Work Step [28]

Description: Click on the column name that is the desired sort parameter

Comments: E.G. Bottle number, state, expiration date, received date, etc.

Work Step [32]

Description: If desired, click on the patient's name or Donor Breast Milk Bottles link and view the details

Facilities:

Facility Name: St. Catharines

Status: Approved

Facility Cerner Owner:

Facility Client Owner: Morosin, Julianne

Authorize Date: Feb 06, 2024

Facility Comments:

Future State: Bridge Transfusion - Bridge Breast Milk - Available Bottles

Cerner Workflow ID: 6974 (v. 10.0) Client Workflow ID: 532

Last updated by Reiser Sanders, Feb 06, 2024 12:24pm (UTC -4 hours)