



Future State: Bridge Transfusion - Bridge Breast Milk - Combine Bottles

Cerner Workflow ID: 6969 (v. 10.0) Client Workflow ID: 534

Last updated by Reiser Sanders, Feb 06, 2024 12:26pm (UTC -4 hours)

Workflow Details:

Workflow Name: Bridge Transfusion - Bridge Breast Milk - Combine Bottles

Workflow State: Future State

Workstream: Other

Venue: Acute Care

Client Owner:

Cerner Owner:

Standard: Yes

Related Workflow(s):

Tags:

Workflow Summary:

Service Line:

Related Solution(s): Bridge Breast Milk Management

Project Name: Niagara Health System:OPT-0297674:NIAG_CD Niagara HIS RFP

TestBuilder Script(s):

Cerner Workflow ID: 6969 (v. 10.0)

Client Workflow ID: 534

Workflow Notes:

Introduced By: WS 3.2

Validated By: WS 5.2

Swim Lane:

Role(s): Nurse (RN)

Department(s):

Security Position(s): Nurse - Neonatology

Women's Health - Nurse

Start/Stop [4]

Description: Access Cerner Bridge

Work Step [6]

Description: Select Combine Bottles

Work Step [8]

Description: Scan the barcode on the bottle label

Decision [68]

Description: Are there additional bottles to be combined?

Future State: Bridge Transfusion - Bridge Breast Milk - Combine Bottles

Cerner Workflow ID: 6969 (v. 10.0) Client Workflow ID: 534

Last updated by Reiser Sanders, Feb 06, 2024 12:26pm (UTC -4 hours)

Decision [17]

Description: Is one bottle being combined into the other bottle?

Comments: If the volume entry field is available and the volumes to be combined from each bottle is different, the bottle with the largest volume to be combined will be the only bottle available to be combined into.

Work Step [30]

Description: Select which bottle is going to have milk combined into it

Comments: Note: It is best practice to combine into the OLDEST bottle

Work Step [23]

Description: If applicable, check the 'Dispose Bottle' button next to all bottles that are getting disposed

Work Step [84]

Description: Click Confirm to update bottle information and to print a new label for the bottle

Start/Stop [87]

Description: Click 'Exit' to return to the patient's chart

Work Step [5]

Description: Scan the patient's wristband or the storage bin label

Decision [16]

Description: Is a new bottle being created?

Work Step [21]

Description: Validate that the 'New Bottle' selection is checked

Work Step [81]

Description: If applicable, check the 'Dispose Bottle' button next to all bottles that are getting disposed

Comments: Can also select "Dispose All" to automatically mark all old bottles for disposal.

Work Step [39]

Description: Click Confirm to print a label for the new bottle

Future State: Bridge Transfusion - Bridge Breast Milk - Combine Bottles

Cerner Workflow ID: 6969 (v. 10.0) Client Workflow ID: 534

Last updated by Reiser Sanders, Feb 06, 2024 12:26pm (UTC -4 hours)

Start/Stop [18]

Description: Click exit and return to the patient's chart

Work Step [76]

Description: Scan the desired printer or select it from the list

Start/Stop [40]

Description: Click 'Exit' to return to the patient's chart

Facilities:

Facility Name: St. Catharines

Status: Approved

Facility Cerner Owner:

Facility Client Owner: Morosin, Julianne

Authorize Date: Feb 06, 2024

Facility Comments: